



FC receptacle

FC receptacle

6,500

12,100

28,700

Ø17

PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF Doric Lenses Inc.. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF Doric Lenses Inc. IS PROHIBITED.

		UNLESS OTHERWISE SPECIFIED:		NAME	DATE	1298, St-Paul, Ancienne-Lorette G2E 1Z3, Canada Tel: 418-877-6886	
		DIMENSIONS ARE IN MM		DRAWN		TITLE:	
		ANGULAR: ± 0.25 deg		CHECKED	M.Vukotic	FRJ_1x1_FC-FC_v04(DRW)	
		ONE PLACE DECIMAL ±0.1		ENG APPR.		REV	
		TWO PLACE DECIMAL ±0.05		MFG APPR.		01	
		THREE PLACE DECIMAL ±0.005		Q.A.		SCALE: 2:1 WEIGHT: 9.64 SHEET 1 OF 1	
		INTERPRET GEOMETRIC TOLERANCING PER:		Third angle projection:		SIZE DWG. NO. REV	
		MATERIAL				B DXXX-YYYY	
NEXT ASSY	USED ON	FINISH					
		BLACK ANODIZED					
APPLICATION		DO NOT SCALE DRAWING					