


# Integrated Fluorescence Mini Cubes (iFMC & iLFMC)



## Differences between generation 1 and 2

Model shown: iFMC6\_IE(400-410)\_E1(460-490)\_F1(500-540)\_E2(555-570)\_F2(580-680)\_S



Model shown: iFMC6-G2\_IE(400-410)\_E1(460-490)\_F1(500-540)\_E2(555-570)\_F2(580-680)\_S

## Improvements

- Amplifiers are now integrated into the cube for a more compact footprint
- Single power input and power switch for multi-detector cubes
- Amplifier controls are located near their corresponding channel
- No more M5 cables