



SEE WEBSITE FOR AVAILABLE LENGTH

LENGTH MAY VARY DEPENDING ON YOUR CONFIGURATION

$\varnothing 0.50$
(1.00mm ALSO AVAILABLE)

4.3 max.

3.2

3.5

ORIENTATION GROOVE

8.0

GRIN

FOCUSING RING

doric

357 rue Franquet, Quebec, QC
G1P 4N7, Canada;
Tel: 418 877 5600;
doriclenses.com;

Dimensions in mm
For Information Purposes Only

PROPRIETARY

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF **Doric Lenses Inc.**. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF **Doric Lenses Inc.** IS PROHIBITED.

DO NOT SCALE DRAWING

TITLE:

TRFC

DWG. NO.

B293-0403-XXX