

Standard	<p><b>Mono fiber</b></p>	<p><b>Two ferrule</b></p>	<p><b>Dual Ferrule</b></p>
	<p>1 fiber / 1 connector</p>	<p>2 fibers / 2 connectors pitch 1.8 to 4.5 mm</p>	<p>2 fibers / 1 connector pitch 0.7 to 1.7 mm</p>
Custom	<p>Array fiber cannula "Bundle-Type" <i>From 2 to 8 fibers</i></p>	<p>Array fiber cannula "Independant" <i>From 2 to 4 fibers</i></p>	<p>High density Fiber optic cannula Array <i>From 9 to 19 fibers</i></p>
	<p>N fibers / 1 connector pitch custom</p>	<p>N fibers / N connectors pitch custom</p>	<p>N fibers / 1 connector pitch hexagonal or custom</p>