



E
400-460 nm ↓

F
500-550 nm ↑

doric
iFMC5

580-650 nm →

SAMPLE ←

IE
400-410 nm ↑

F

BNC CONNECTOR

POWER SUPPLY JACK 2.5 MM

ON/OFF SWITCH

AMPLIFICATION
MODE SWITCH

GAIN LEVEL
SWITCH

FC RECEPTACLE

Dimensions in mm
For Information Purposes Only

PROPRIETARY
THE INFORMATION CONTAINED IN THIS DRAWING
IS THE SOLE PROPERTY OF Doric Lenses Inc..
ANY REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
Doric Lenses Inc. IS PROHIBITED.

DO NOT SCALE DRAWING

doric 357 rue Franquet, Quebec, QC
G1P 4N7, Canada;
Tel: 418 877 5600;
doriclenses.com;

TITLE:
iFMC5 1DET 3FC

DWG. NO.
B340-1254-001

4-6-2020