



doric

i1FMC6

IE 400-410nm E1 460-490nm E2 555-570nm F1 500-540nm F2 580-680nm

SAMPLE

BNC CONNECTOR

GAIN LEVEL SWITCH

ON/OFF SWITCH

FC RECEPTACLE

POWER SUPPLY JACK 3.5 MM

doric

357 rue Franquet, Quebec, QC
G1P 4N7, Canada;
Tel: 418 877 5600;
doriclenses.com;

TITLE:

**i1FMC6-G3 2 détecteurs
- 3 DELs**

DWG. NO.

B340-1332-XXX

Dimensions in mm
For Information Purposes Only

PROPRIETARY
THE INFORMATION CONTAINED IN THIS DRAWING
IS THE SOLE PROPERTY OF Doric Lenses Inc..
ANY REPRODUCTION IN PART OR AS A WHOLE
WITHOUT THE WRITTEN PERMISSION OF
Doric Lenses Inc. IS PROHIBITED.

DO NOT SCALE DRAWING

8-4-2022